

## Polarizers

# Opto-Physics Infrared Polarizers (GS57500 Series)

Polarizers are commonly used to polarize radiation from unpolarized sources, attenuated radiation from polarized sources, or act as polarizing beamsplitters. Specac offers a range of holographic wire polarizers laid onto a transmitting substrate material for use in the 2 - 35 $\mu\text{m}$  (5000 $\text{cm}^{-1}$  - 285 $\text{cm}^{-1}$ ) spectral range.

These precision polarizers are manufactured in a class 1000 clean room facility at Specac's United Kingdom factory, by means of a holographic fabrication technique originally developed in conjunction with the United Kingdom's National Physical Laboratory (NPL).

The process involves exposing a photo-resist coating on a suitable material substrate to an interferometrically-generated fringe pattern from a monochromatic UV source. The regular sinusoidal profile of the developed photo-resist is subsequently metal coated at an oblique angle to create an array of fine parallel lines at a set period.

This technique lends itself well to the generation of extremely uniform sub-micron grid wire spacings at 2500 lines / mm, which have significantly reduced level of light scattering in comparison to traditional ruled wire grid polarizers. As the wire grid is formed on the photo-resist itself, the technique is also well suited to fabricating polarizers on substrates that do not otherwise lend themselves to the ruling process.

Specac offers a wide range of polarizers on infrared material substrates such as Barium Fluoride (BaF<sub>2</sub>), Calcium Fluoride (CaF<sub>2</sub>), KRS-5, Zinc Selenide (ZnSe), and Germanium (Ge), in a range of categories to meet a broad scope of customer requirements.

Illustrations, descriptions and specifications in this datasheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time.

For the latest details please contact your local Specac office or representative.



### Features of GS57500 Series Polarizers

- 2500 lines/mm on substrate
- Choice of KRS-5, CaF<sub>2</sub>, BaF<sub>2</sub>, and ZnSe substrates
- 18mm C.A. / 25mm O.D. and 34mm C.A. / 50mm O.D. options
- Free standing by use of own ring mount
- Economical price


### Applications

- Infrared spectroscopy of materials (typically plastics / polymers and crystallography)
- Infrared microscopy (sample characterization)
- NIR/Mid-IR thermal imaging systems
- Plasma diagnostics
- Beamsplitters in polarized light interferometry
- Analysis in infrared astronomy
- Low power laser polarization and beam attenuation
- Coupling devices for Mid-IR and long wavelength lasers

## Polarizers

# Opto-Physics Infrared Polarizers (GS57500 Series)

Typical Transmission (K1) and Attenuation (K2) Plots for Opto-Physics Polarizers



Substrate Material	Wavelength (μm) / Wavenumber (cm <sup>-1</sup> )	Typical Extinction Ratio (K1/K2)
KRS-5	2.5 / 4000	100
	5 / 2000	300
	10 / 1000	300
CaF2	2.5 / 4000	100
	5 / 2000	300
	10 / 1000	300
BaF2	2.5 / 4000	100
	5 / 2000	300
	10 / 1000	400
ZnSe	2.5 / 4000	100
	5 / 2000	300
	10 / 1000	400

### Opto-Physics Polarizer Specifications

Wire Grid Spacing (lines/mm)	2500	2500
Mounting Ring Diameters (mm)	25.0 +0.0/-0.2	50.0 +0.0/-0.2
Mounting Ring Thickness (mm)	5.0 ±0.1	6.0 ±0.1
Polarizer Clear Aperture (mm)	18.0 ±0.1	34.0 ±0.1
Substrate Thickness (mm)	2.0 +0.0/-0.2	4.0 +0.0/-0.2

	Clear Aperture (mm)	Mounting Ring Dia. (mm)	Parallelism (arc minutes)	Flatness (fringes per inch @633 nm)	Part Number
KRS-5	18	25	5	4	GS57500
	34	50	5	4	GS57504
CaF2	18	25	3	2	GS57501
	34	50	3	2	GS57505
BaF2	18	25	3	2	GS57502
	34	50	3	2	GS57506
ZnSe	18	25	3	2	GS57503
	34	50	3	2	GS57507

SPECAC LTD.,  
 River House 97 Cray Avenue, Orpington,  
 Kent. BR5 4HE UK  
 T: +44 (0) 1689 873134 F: +44 (0) 1689 878527  
 E: sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.,  
 301 Berkeley Drive Suite B, Swedesboro, NJ 08085 USA  
 T: +1 856 241 1925 F: +1 856 241 1926  
 E: sales@specac.com W: www.specac.com  
 Specac is part of Smiths Group plc

MD 15 - 05/13